




SONICA® 3200

6 litres

Ultrasonic cleaning tank available in version M, MH, ETH and EP.
Wide tank suitable for cleaning dental instruments, watch/clock parts, small printed circuit boards, hybrid circuits and precision mechanical parts

Tank size

300x240x100mm



General informations

- Average mains power when not using heater 180W
- Average mains power when using heater 355W
- HF Power 720W p-p.
- Ultrasound frequency 45Khz \pm 5%

- Standard drain
- Number of transducers 4
- Capacity 6lts.
- Case and tank made in stainless steel AISI 304

Control panel



VERSION M
The simplest version in the SONICA® range is equipped with a timer which can be set for to 15 minutes. It can also be set to manual operation.



VERSION ETH
Cleaning times and temperatures controlled electronically by a microprocessor.
Waterproof control panel.



VERSION MH
This is equipped with a timer which can be set for up to 15 minutes and can also be set to manual operation. The thermostat for the heater is set at 60°C.



VERSION EP
The most technologically advanced of the SONICA® range.

- New Sweep system generator.
- Modern digital display control panel windows.
- Degas time setting.

- Cleaning times from 0 to 99 minutes.
- Tank liquid temperature setting up to 70°C.
- Memory for 9 cleaning programmes.
- Jet program function to start quickly and easily cleaning cycle.

Accessories

- Stainless steel holder basket
- Small stainless steel immersion basket
- Glass beaker
- Plastic beaker
- Plastic positioning lid for 4 beaker
- Stainless steel lid
- Stainless steel cassettes carrier
- SONICA® stainless cassettes carrier

Liquid detergents

SONICA® UG Concentrated detergent for general use. Bottle Kg1.

SONICA® SF Concentrated detergent for removing heavy grease from very dirty spots. Bottle Kg1.

SONICA® AF Concentrated detergent for removing limestone deposit stains and large quantities of oxides and rust. Bottle Kg1.

SONICA® AC Concentrated detergent for removing limestone traces and small amounts of oxides and rust. Bottle Kg1.

SONICA® RC Ready-to-use detergent for removing cement deposits from dental instruments. Bottle Kg1.

SONICA® RA Ready-to-use detergent for removing alginate deposits from impression trays in dental application. Bottle Kg1.

SONICA® RG Special neutral detergent for removing plaster deposits. Bottle Kg1.

