



SONICA® 4300

18 litres

Medium-sized ultrasonic cleaning tank available in version M, MH, ETH and EP. Suitable for cleaning medium-sized objects such as small moulds, cleaning sieves, filters and endoscopes. This model can accommodate three large dental instrument trays with stainless steel holders.

Tank size

330x300x200mm



General informations

- Average mains power when not using heater 300W
- Average mains power when using heater 800W
- HF Power 1200W p-p.
- Ultrasound frequency 40Khz \pm 5%
- Standard drain and ball valve
- Number of transducers 4
- Capacity 18lts.
- Case and tank made in stainless steel AISI 304

Control panel



VERSION M
The simplest version in the SONICA® range is equipped with a timer which can be set for to 15 minutes. It can also be set to manual operation.



VERSION ETH
Cleaning times and temperatures controlled electronically by a microprocessor.
Waterproof control panel.



VERSION MH
This is equipped with a timer which can be set for up to 15 minutes and can also be set to manual operation. The thermostat for the heater is set at 60°C.



VERSION EP
The most technologically advanced of the SONICA® range.

- New Sweep system generator.
- Modern digital display control panel windows.
- Degas time setting.
- Cleaning times from 0 to 99 minutes.
- Tank liquid temperature setting up to 70°C.
- Memory for 9 cleaning programmes.
- Jet program function to start quickly and easily cleaning cycle.

Accessories

- Stainless steel holder basket
- Small stainless steel immersion basket
- Glass beaker
- Plastic beaker
- Plastic positioning lid for 4 beaker
- Stainless steel lid
- Stainless steel cassettes carrier
- Special stainless steel holder basket

Liquid detergents

SONICA® UG Concentated detergent for general use. Bottle Kg1.

SONICA® SF Concentated detergent for removing heavy grease from very dirty spots. Bottle Kg1.

SONICA® AF Concentated detergent for removing limestone deposit stains and large quantities of oxides and rust. Bottle Kg1.

SONICA® AC Concentated detergent for removing limestone traces and small amounts of oxides and rust. Bottle Kg1.

SONICA® RC Ready-to-use detergent for removing cement deposits from dental instruments. Bottle Kg1.

SONICA® RA Ready-to-use detergent for removing alginate deposits from impression trays in dental application. Bottle Kg1.

SONICA® RG Special neutral detergent for removing plaster deposits. Bottle Kg1.

